

# VITRO

blood tube



## Use case



### **Blood tube production with plastic as the glass substitute material**

ZAHORANSKY also relies 100 percent on plastic as a glass substitute in the manufacture of blood collection tubes. This is precisely where its advantages – no risk of breakage, easier handling, uncomplicated storage – come into their own. ZAHORANSKY meets all requirements for high-precision manufactured products with many years of experience in mold manufacture. Thanks to sites distributed worldwide, the corresponding production capacities can be provided flexibly, quickly and, above all, close to the customer. This not only benefits a company's ecological footprint, but thanks to short distances, it also ensures faster and cheaper on-site service calls and the possibility of direct support.

ZAHORANSKY currently has more than 18 injection molding machines available. Each tool in it is designed for 100 million parts per year in different sizes. The systems can be used for trial installation, among other things. In the field of subsequent automation, ZAHORANSKY works with a partner who completes the blood collection tubes. Thanks to the "one-stop shopping" concept, customers receive everything from a single source.